

SeeMore™ USB Probes



We revolutionize imaging devices.

VETERINARY PROBES



GV 5.0 MHz

- Center Frequency: 5.0 MHz
- Pulse Frequency: 5.0 MHz and 3.5 MHz
- Scan Angle: 60°
- Focal Point: 6 cm
- Scan Depth (Max): 15 cm
- Applications : Large Animals; Abdomen (kidney, liver, gall bladder), Uterus, Bladder

SV 7.5 MHz

- Center Frequency: 7.5 MHz
- Pulse Frequency: 5.0 MHz to 24 MHz
- Scan Angle: 60°
- Focal Point: 2.0 cm
- Scan Depth (Max): 10 cm
- Applications : Small Animals; Abdomen (kidney, liver, gall bladder), Uterus, Bladder

MINIMUM SYSTEM REQUIREMENTS

- Computer Operating System: Windows XP, Vista, or Mac (running Windows)
- Processor: 1 GHz Minimum
- Memory: 512 K (recommended 1 GB or more)
- USB 2.0 port
- Minimum Display: 1024 x 600 with 32 bit color