

SeeMore™ USB Probes



We revolutionize imaging devices.

ENDOCAVITY PROBES



EC 7.5 MHz

EC 7.5 MHz

- Center Frequency: 7.5 MHz
- Pulse Frequency: 5.0 MHz to 12 MHz
- Scan Angle: 90°
- Focal Point: 2.5 cm
- Scan Depth (Max): 10 cm
- Clinical Applications :
Transrectal: Prostate
Transvaginal: Uterus, Adnexa
- Notes: End Fire, Optional Biopsy Guide, Prostate Volume Calculation